



Bioinformatics Summer School

Long-reads Transcriptomics

Ana Conesa & Fabian Jetzinger

I²SysBio, Valencia, Spain
BioBam Bioinformatics, Valencia, Spain

- 1 Biological replicates and transcriptome reconstruction
- 2 Differential Expression vs. Differential Isoform Usage
- 3 tappAS

Section 1

Biological replicates and transcriptome reconstruction

Understanding the biological replicates
influence transcriptome reconstruction

- Distinguish general patterns from sample-specific noise
- Detect and discard outliers, low-quality samples
- Increase power of statistical analyses
- Enable a wider range of tools


Join & Call

- Combine all reads → transcriptome reconstruction
- Detect lowly-expressed transcript variants




Call & Join

- Transcriptome reconstruction on each sample → combine annotations (e.g. TAMA Merge)
- Detect sample-specific transcript variants



Results of both strategies vary across transcriptome reconstruction tools!

Transcriptome diversity




Structural Category



Transcriptome diversity varies based on strategy and reconstruction tool


Transcriptome
diversity




Structural Category



Assignment of
reads




Transcriptome diversity



Structural Category



Assignment of reads



Takeaways:

Performance of strategy depends on reconstruction tool

Majority of reads assigned to FSM (reference, known) transcripts

Consider choice of strategy and tool carefully based on research objectives


Section 2

Differential Expression vs. Differential Isoform Usage

Understanding differential expression at isoform resolution

Differential Expression


- Gene- or transcript-level
- Compares **total expression** across conditions




DGE

Differential Expression

- Gene- or transcript-level
- Compares **total expression** across conditions




DGE




DTE

Differential Expression


- Gene- or transcript-level
- Compares **total expression** across conditions



DGE



DTE




DTE

Differential Expression

- Gene- or transcript-level
- Compares **total expression** across conditions


Differential Isoform Usage

- Compares **relative expression of a gene's isoforms** across conditions



Differential Isoform Usage

- Compares **relative expression of a gene's isoforms** across conditions




DGE



DIU (major isoform switch)

Differential Isoform Usage


- Compares **relative expression of a gene's isoforms** across conditions



DGE



DIU (major isoform switch)



DGE



DIU



- Other, similar types of analyses exist, e.g. Differential Exon Usage
- Popular tools: edgeR, DESeq2, DEXSeq, NOISeq (can simulate replicates), DRIMSeq, etc.
- Traditional tools (edgeR, DESeq2, etc.) have been shown to also work well in long reads

See:

Dong, Xueyi et al. "Benchmarking long-read RNA-sequencing analysis tools using in silico mixtures." *Nature methods* vol. 20,11 (2023): 1810-1821. doi:10.1038/s41592-023-02026-3

How would you describe the following situation of Gene A and its 3 isoforms?

- Differential Gene Expression?
- Differential Transcript Expression?
- Differential Transcript Usage?
- Major Isoform Switch?



	Tr. A.1	Tr. A.2	Tr. A.3
Condition 1	1000	600	400
Condition 2	500	100	400


Section 3

tappAS

Applying differential expression and usage analyses
and much more...

- Gene-, transcript- or CDS-level differential expression
- Transcript- or CDS-level differential usage
- Gene Set Enrichment Analysis (GSEA)
- Visualization of functional annotation

tappAS project configuration



New Project

Project Information

Name: isotools_isoannotlite

Annotation Data

Species: Homo sapiens

☐ Use application's annotation features file - requires download first time used

Select annotation features file

☒ Specify annotation features file's location

/home/train/longTREC/day2/data/sqanti3/data/h1_endo_chr8_isotools_isoannot.gff3

Experiment Design and Data

Experiment Type: Two-Group Comparison

Design File: /home/train/longTREC/day2/data/sqanti3/data/exp_design.txt

Matrix Data File: /home/train/longTREC/day2/data/sqanti3/data/h1_endo_chr8_all_count.txt

Note: First group in design file must be the control and expression matrix MUST contain raw counts - see Help

☒ Normalize data by TMM method (recommended for raw counts)

☒ Expression Values Filter (recommended):

Low count values cutoff: 1.0 CPM Coefficient of variation cutoff: 100.0 %

Transcripts Filter: Do not filter transcripts

Use button on right or type in full file path and name

Help Cancel OK

Thank You!



For more information about the LongTREC Summer School:

<https://longtrec.eu>